

■表 8 Multiplex PCR による pAmpC β-ラクタマーゼ遺伝子検索用プライマー

pAmpC β-ラクタマーゼ遺伝子	プライマー配列	PCR 増幅産物サイズ (bp)
MOX type	5'-GCTGCTCAAGGAGCACAGGAT-3' 5'-CACATTGACATAGGTGTGGTGC-3'	520
CIT type	5'-TGGCCAGAACTGACAGGCAAA-3' 5'-TTTCTCCTGAACGTGGCTGGC-3'	462
DHA type	5'-AACTTTCACAGGTGTGCTGGGT-3' 5'-CCGTACGCATACTGGCTTTGC-3'	405
ACC type	5'-AACAGCCTCAGCAGCCGGTTA-3' 5'-TTCGCCGCAATCATCCCTAGC-3'	346
EBC type	5'-TCGGTAAAGCCGATGTTGCGG-3' 5'-CTTCCACTGCGGCTGCCAGTT-3'	302
FOX type	5'-AACATGGGGTATCAGGGAGATG-3' 5'-CAAAGCGCGTAACCGGATTGG-3'	190