

■表 11 主な菌種における PFGE 解析に適切な制限酵素の例<sup>5)</sup>

被検菌	制限酵素
<i>E. coli</i>	XbaI, BlnI (AvrII), NotI, SfiI
<i>Klebsiella</i> spp.	XbaI, AsnI
<i>Enterobacter</i> spp.	XbaI, SpeI
<i>P. mirabilis</i>	SfiI, NotI
<i>Citrobacter</i> spp.	XbaI
<i>S. marcescens</i>	XbaI
<i>Salmonella</i> spp.	XbaI, BlnI (AvrII), NotI
<i>Shigella</i> spp.	XbaI, BlnI (AvrII), SfiI
<i>P. aeruginosa</i>	SpeI, XbaI, DraI, SspI
<i>Acinetobacter</i> spp.	SmaI, ApaI
<i>H. influenzae</i>	SmaI, RsrII
<i>Staphylococcus</i> spp.	SmaI, CspI, SstII, SgrAI
<i>Enterococcus</i> spp.	SmaI, ApaI
<i>S. pneumoniae</i>	SmaI, ApaI
<i>Streptococcus</i> spp. (group A and B)	SmaI
<i>L. pneumophila</i>	SfiI, NotI
<i>C. difficile</i>	SmaI, SacII, SstII, NruI
<i>Camphylobacter</i> spp.	SmaI, SalI, KpnI